



Absorbency: 6.7 oz/198 mL, size Medium



## Booster Pads

A soft, flow-through pad designed to fill to capacity and then pass additional fluid into the primary garment.

Product Features:

BR

KG

CL

MH

WI

RT


AS

LF

ITEM #	SIZE	LENGTH/WIDTH	CAPACITY	PACKAGING	HCPC
2770	Small	12" x 3.25"	5.5 oz/163 mL	6 bags of 30 (180)	T4535
2760	Medium	12" x 4.25"	6.7 oz/198 mL	8 bags of 25 (200)	T4535
2762	Extended	15" x 4.25"	8.0 oz/236 mL	8 bags of 25 (200)	T4535


## What is a booster pad?

Booster pads have a flow-through design that allows the pad to fill to capacity first and then pass additional fluid to the host (primary) garment. Booster pads are placed inside any disposable undergarment with a moisture barrier, boosting absorbency to extend wear-time. They are also easy to change and dispose.




## Alternative Placements for Booster Pads

- Position booster pad horizontally across the stomach, adhesive side up. Apply absorbent undergarment as usual and adhere pad.
- Strategically place the booster pad over the hip to help prevent leaking for side sleepers.
- For women, position the pad toward the back. For men, position the pad toward the front.
- For men, DO NOT remove adhesive backing. Fold booster pad over top of penis with release paper to the outside. Apply absorbent undergarment as usual.



Tummy pad placement



Fold pad in half