

SEAT OPTIONS FOR THE EZEE LIFE REHAB SHOWER COMMODES

EZee Life™ now offers up to three different seat options for its aluminum rehab shower commodes, depending on the model, to meet the needs of a large range of users and caregivers.

Our commodes are fitted with a carefully chosen premium quality padded seat, back and arms to enhance client comfort and postural support while being able to withstand the rigours of use and vigorous cleaning solvents found in institutions and homes.

All our seats have a unique tool-less design that can be easily removed from the frame for cleaning while still providing durability and stability. The ease of cleaning these rehab commodes make them one of the most effective commodes available for situations where infection control is paramount.

OPEN FRONT SEAT

The open front commode seat is the standard seat found on all EZee Life™ aluminum rehab shower commodes. This seat allows hygiene access by the user or a caregiver from the front of the commode and has a padded insert that fills in the commode opening for those who don't want or need to use the commode feature.

The open front seat is available on all models and sizes.



CLOSED FRONT SEAT


The closed front commode seat is an optional seat that fits several models of EZee Life™ aluminum rehab shower commodes. The closed front is for users who do not require hygiene access while seated. This seat comes with a full padded top cover for users who don't need to have a commode opening all the time but want the option to switch back and forth.

The closed front seat is available on 180, 185, 190 and 195.



4-WAY SEAT

The 4-Way Seat is a standard open front commode seat that can be mounted onto the commodes facing the front, rear, left or right side. This seat option allows hygiene access from any direction. The seat can be changed from client to client and from caregiver to caregiver depending on preference and needs.



The 4-Way Seat is only available on the 180 and 190 commode models.



4-WAY FRONT



4-WAY LEFT



4-WAY REAR



4-WAY RIGHT