

# ProBasics® SatinAir™ Bariatric Mattresses

HCPCS CODE: E0277

The SatinAir Bariatric Mattress Series is the first choice of clinicians to treat and prevent bariatric users from pressure ulcers.

- Alternating pressure cycles of 10, 15, 20 & 25 allows clinicians to customize pressure redistribution
- Low air loss function controls and manages moisture to prevent skin maceration
- Cell-on-cell mattress design prevents user from bottoming out in case of power failure
- 10" air cells provide comfort and support
- Audible low pressure alarms alert clinicians when preset pressure levels decrease
- Static mode allows for easy transfer in and out of bed
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode
- CPR pull tab quickly releases air from mattress



Model #	Overall Dimensions	Weight Capacity	Compressor	Alternating Time Cycle	Alarms	Regulation	Limited Warranty
PBAPMBARI42	80"L 42"W 10"H	650lb.	15 LPM	10, 15, 20, 25 Minute Intervals	Audio	UL Approved	2 Years
PBAPMBARI48	80"L 48"W 10"H	650lb.	15 LPM	10, 15, 20, 25 Minute Intervals	Audio	UL Approved	2 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

# C SUPPORT