

All-in-one brief

Abri-Form Premium



Product features

- Breathable brief to maintain healthy skin and optimum comfort for the user.
- Top Dry-system, ensuring daily comfort and skin integrity.
- High quality flexible tapes for easy and optimal fit and easy repositioning.
- Barriers to ensure the user protection against leakage.
- Odour control for discretion.
- Wetness indicator with graduation scale.



Application guide



Fold the brief lengthways with backsheet turned outwards. Do not touch the inside of the brief.



Place the wearer on their side. Pull the brief from front to back.



Place the wearer on their back. Fan out the front of the brief, creating a channel shape between the legs.



Fix bottom tapes on both sides for a snug, comfortable fit.



Fix top tapes angled downwards over the pleat.



Abri-Form Junior

Abri-Form Junior is designed for juniors with quick and uncontrolled urination, and who would otherwise fall between baby diapers and adult products. With all premium quality features and security.



Abri-Form Premium

Incontinence

Bed Protection

Protective Wear

Wound Care

Body Wash

Skin Care

Gloves

■ Abri-Form Junior, Premium

Art. no.	Product/ Colour Code	Hip size (cm)	Absorption Rothwell	Drops	Units per pack/ ctn
43050	Junior XS2	40-60	1500 ml	7	32/128










■ Abri-Form, Premium

Art. no.	Product/ Colour Code	Hip size (cm)	Absorption Rothwell	Drops	Units per pack/ ctn
43054	XS2	50-60	1400 ml	7	32/128
43055	S2	60-85	1800 ml	7	28/84
43056	S4	60-85	2200 ml	9	22/66










■ Abri-Form, Premium

Art. no.	Product/ Colour Code	Hip size (cm)	Absorption Rothwell	Drops	Units per pack/ ctn
43049	M0	70-110	1500 ml	5	26/104
43061	M1	70-110	2000 ml	6	26/104
43060	M2	70-110	2600 ml	7	24/96
43062	M3	70-110	3100 ml	8	22/88
43063	M4	70-110	3600 ml	9	14/56










■ Abri-Form, Premium

Art. no.	Product/ Colour Code	Hip size (cm)	Absorption Rothwell	Drops	Units per pack/ ctn
43059	L0	100-150	2000 ml	5	26/104
43066	L1	100-150	2500 ml	6	26/104
43065	L2	100-150	3100 ml	7	22/88
43067	L3	100-150	3400 ml	8	20/80
43068	L4	100-150	4000 ml	9	12/48










Abri-Form Premium

■ Abri-Form, Premium

	Art. no.	Product/ Colour Code	Hip size (cm)	Absorption Rothwell	Drops	Units per pack/ ctn
	43069	XL2	110-170	3400 ml	7	20/80
	43071	XL4	110-170	4000 ml	9	12/48

ABRI • FORM



■ Abri-Form, Premium

	Art. no.	Product/ Colour Code	Hip size (cm)	Absorption Rothwell	Drops	Units per pack/ ctn
	4730	M1	70-110	2000 ml	6	10/60
	4740	M2	70-110	2600 ml	7	10/60
	4735	L1	100-150	2500 ml	6	10/60
	4745	L2	100-150	3100 ml	7	10/60

ABRI • FORM

