

Underpads



Simplicity™ Extra Underpad



- Spunbond non-woven topsheet allows fluid to pass quickly into the absorbent core
- Polypropylene backsheet protects against leakage and resists melting

Code #	Size	Absorbent Core	Backsheet Color	Packaging	Ship Case
949	17" x 24"	Fluff	Lt. Blue	Bulk	300
1038	23" x 24"	Fluff	Lt. Blue	50/Bag, 4 Bags/Case	200
6349	23" x 24"	Fluff	Lt. Blue	10/Bag, 20 Bags/Case	200
1093	23" x 36"	Fluff	Lt. Blue	10/Bag, 15 Bags/Case	150
949B10	30" x 30"	Fluff	Lt. Blue	10/Bag, 10 Bags/Case	100



Simplicity™ Basic Underpad



- Spunbound non-woven topsheet allows fluid to pass quickly into the absorbent core
- Polypropylene backsheet protects against leakage and resists melting

Code #	Size	Absorbent Core	Backsheet Color	Packaging	Ship Case
7105	17" x 24"	Fluff	Lt. Blue	Bulk	300
7134	23" x 24"	Fluff	Lt. Blue	50/Bag, 4 Bags/Case	200
7136	23" x 24"	Fluff	Lt. Blue	10/Bag, 20 Bags/Case	200
7174	23" x 36"	Fluff	Lt. Blue	50/Bag, 3 Bags/Case	150
7176	23" x 36"	Fluff	Lt. Blue	10/Bag, 15 Bags/Case	150

Specialty Underpads



Wings™ Quilted Premium Breathable Underpad



- Strong cloth-like backsheet resists tearing while maintaining the pad's integrity
- Breathable backsheet allows air to freely circulate and allows patients to garner the benefits from an air therapy bed
- Super absorbent core is designed to quickly lock in fluid and neutralize unpleasant odors
- Moisture vapor permeable film allows moisture vapor (gas) to transfer away from the patient, while protecting bed linens/clothing from bodily fluids, leaving the patient drier and more comfortable

Code #	Size	Absorbent Core	Backsheet Color	Packaging	Ship Case
P2336B	23" x 36"	Airlaid	Light Blue	12/bag, 6 bags/case	72
P3030B	30" x 30"	Airlaid	Light Blue	10/bag, 6 bags/case	60
P3036B	30" x 36"	Airlaid	Light Blue	10/bag, 4 bags/case	40



Wings™ Breathable Plus Underpad



- Breathable backsheet allows air to freely circulate and allows patients to garner the benefits of a low air loss bed
- Contains super-absorbent polymer for maximum fluid retention while reducing odor and keeping wetness away from skin

Code #	Size	Absorbent Core	Backsheet Color	Packaging	Ship Case
988	23" x 36"	Fluff/Polymer	Lt. Blue	12/Bag, 6 Bags/Case	72
988B6	23" x 36"	Fluff/Polymer	Lt. Blue	6/Bag, 12 Bags/Case	72
988HD	23" x 36"	Fluff/Polymer	Lt. Blue	12/Bag, 6 Bags/Case	72
984	30" x 36"	Fluff/Polymer	Lt. Blue	10/Bag, 6 Bags/Case	60



Wings™ Specialty Underpad



- Extra-large size provides maximum coverage, absorption and protection equivalent to three or four regular fluff underpads
- Contains super-absorbent polymer for maximum fluid retention while reducing odor and keeping wetness away from skin
- Adhesive strips/tuckable wings minimize pad movement and bunching

Code #	Size	Absorbent Core	Backsheet Color	Packaging	Ship Case
Adhesive Underpad					
959	30" x 36"	Fluff/Polymer	Lt. Blue	12/Bag, 6 Bags/Case	72
Tuckable Underpad					
995	36" x 70"	Fluff/Polymer	Peach	12/Bag, 4 Bags/Case	48