

Bariatric Sentra Extra-Heavy-Duty Wheelchair

Dual Cross Brace



HCPCS
K0007

- Weight Capacity: 500 lbs.
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable, heavy-gauge reinforced naugahyde upholstery.
- Comes standard with carry pocket on backrest.
- Reinforced side frame, caster journals and dual cross bars.
- Dual-position front fork.
- Dual-position front frame.
- Composite Mag-style wheels with chrome hand rim are lightweight and maintenance free.
- Hemi-adjustable caster is optional (Part #STDS1074).
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Padded armrests provide added patient comfort.
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings.
- Comes with push-to-lock wheel locks.

Item #	Seat Width	Armrest	Product Front Riggings	UOM
STD20DDA-SF	20"	Detachable Desk Arm	Swing-away Footrests	1/cs
STD20DDA-ELR	20"	Detachable Desk Arm	Swing-away, Elevating Leg rests	1/cs
STD20DFA-SF	20"	Detachable Full Arm	Swing-away Footrests	1/cs
STD20DFA-ELR	20"	Detachable Full Arm	Swing-away, Elevating Leg rests	1/cs
STD20ADDA-SF	20"	Detachable and Adjustable Height Desk Arm	Swing-away Footrests	1/cs
STD20ADDA-ELR	20"	Detachable and Adjustable Height Desk Arm	Swing-away, Elevating Leg rests	1/cs
STD20ADFA-SF	20"	Detachable and Adjustable Height Full Arm	Swing-away Footrests	1/cs
STD20ADFA-ELR	20"	Detachable and Adjustable Height Full Arm	Swing-away, Elevating Leg rests	1/cs
STD22DDA-SF	22"	Detachable Desk Arm	Swing-away Footrests	1/cs
STD22DDA-ELR	22"	Detachable Desk Arm	Swing-away, Elevating Leg rests	1/cs
STD22DFA-SF	22"	Detachable Full Arm	Swing-away Footrests	1/cs
STD22DFA-ELR	22"	Detachable Full Arm	Swing-away, Elevating Leg rests	1/cs
STD22ADDA-SF	22"	Detachable and Adjustable Height Desk Arm	Swing-away Footrests	1/cs
STD22ADDA-ELR	22"	Detachable and Adjustable Height Desk Arm	Swing-away, Elevating Leg rests	1/cs
STD22ADFA-SF	22"	Detachable and Adjustable Height Full Arm	Swing-away Footrests	1/cs
STD22ADFA-ELR	22"	Detachable and Adjustable Height Full Arm	Swing-away, Elevating Leg rests	1/cs
STD24DDA-SF	24"	Detachable Desk Arm	Swing-away Footrests	1/cs
STD24DDA-ELR	24"	Detachable Desk Arm	Swing-away, Elevating Leg rests	1/cs
STD24DFA-SF	24"	Detachable Full Arm	Swing-away Footrest	1/cs
STD24DFA-ELR	24"	Detachable Full Arm	Swing-away, Elevating Leg rests	1/cs
STD24ADDA-SF	24"	Detachable and Adjustable Height Desk Arm	Swing-away Footrests	1/cs
STD24ADDA-ELR	24"	Detachable and Adjustable Height Desk Arm	Swing-away, Elevating Leg rests	1/cs
STD24ADFA-SF	24"	Detachable and Adjustable Height Full Arm	Swing-away Footrests	1/cs
STD24ADFA-ELR	24"	Detachable and Adjustable Height Full Arm	Swing-away, Elevating Leg rests	1/cs

specifications

Width Open with 20" Seat	29"	Product Weight with 20" Seat	54 lbs.
Width Open with 22" Seat	31"	Product Weight with 22" Seat	56 lbs.
Width Open with 24" Seat	32.5"	Product Weight with 24" Seat	66.25 lbs.
Width Closed	13"	Carton Shipping Weight with 20" Seat	60 lbs.
Height	37"	Carton Shipping Weight with 22" Seat	62 lbs.
Seat-To-Floor Height	17.5" - 19.5"	Carton Shipping Weight with 24" Seat	72 lbs.
Back Height	18"	Weight Limit	500 lbs.
Seat Depth	18"	Warranty	Limited Lifetime
Length with Riggings	42"		

accessories

Telescoping IV Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

Articulating Elevating Legrests

STDELR-AL	Articulating Legrest For Manual Wheelchairs	1/cv
-----------	---	------

Universal Forward Stabilizer

STDS860	Universal Forward Stabilizer	1/pr
STDS862	Universal Forward Stabilizer	1/pr

Anti Tippers

STDS802	Anti-Tippers With Wheels	1/pr
---------	--------------------------	------

6" Wheel Lock Extension

STDS801	Wheel Lock Extension	1/ea
---------	----------------------	------

Seat Rail Extension

STDS2420	Seat Rail Extension Kit, 24" models only	1/pr
STDS811	Seat Rail Extensions	1/pr

Bariatric Seat Belts

STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric	1/ea

Universal Oxygen "E" Cylinder Holder

STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea

Overhead Anti-Theft Device

STDS823	Overhead Anti-Theft Device	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea

Anti-Fold Device

STDS806	Anti-Fold Bar	20" and 22" Only	1/ea
---------	---------------	------------------	------

Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

Heel Loops

STDS831	Heel Loops	1/bx
---------	------------	------

Limb Support

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea

Chart Carry Pocket

STDS836	Chart Carry Pocket	20" Wheelchairs	1/ea
---------	--------------------	-----------------	------