



KST001



A




Flame Red



Flame Blue

10289RD



A

B

C

D

K. D. Handy Trolley, 4" Casters

Deluxe 3 Wheel Aluminum Rollator, 7.5" Casters

Loop Locks

Item #	UOM
KST001	1/cs

Item #	Color	UOM
10289RD	Flame Red	1/cs
10289BL	Flame Blue	1/cs

- Angled handles provide for an improved grip and stability.
- Comes with 2 removable trays.
- Tool-free removable casters with euro-style clips. (Figure A)
- The bottom tray is set forward making it easy and safe to walk with.
- Durable steel construction.

- Comes standard with basket, tray and pouch. (Figure A)
- Basket conveniently folds for easy storage. (Figure B)
- Ships in ultra compact carton. Comes with frame folded and compact. Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob. (Figure C)
- 7.5" casters for indoor or outdoor use.
- Handles easily adjust in height with self-threading knob. (Figure C)
- Special loop lock made of internal aluminum casting operates easily and ensures safety. (Figure D)
- Easy, one-hand folding.

specifications

Table/Tray Dimensions	18" (W) x 19" (D)
Base Dimensions	18" (W) x 19" (D)
Height	27.5" - 32.5"
Weight	12 lbs.
Carton Shipping Weight	13 lbs.
Warranty	Limited Lifetime

specifications

Handle Height	31" - 38"
Unit Dimensions	25" (W) x 24" (L)
Carton Shipping Weight	17 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

accessories

KST001-TY	Replacement Tray for KST001 and KST001HB
KST001C	Caster

accessories

815B	Replacement Basket	1/ea
------	--------------------	------