



## All-Terrain Knee Walker

Item #	Color	UOM
990X	Silver Vein	1/ea

- Two 16" front wheels and two 12" rear wheels with enhanced traction effortlessly maneuver over gravel, dirt, grass, mulch, and uneven pavement
- Heavy-duty, extended-wear disc brake enables quick stopping for unparalleled safety
- Knee cushion with suspension is crafted for comfort to reduce potential pressure on shin and knee
- Foldable steering column allows for secure, compact storage and transport
- Complimentary water bottle and custom bottle holder make staying hydrated during travel easy

specifications	
Handle Height	32.875" – 37.625"
Knee Pad Height	18.375" – 23.25"
Pad Dimensions	14" x 6.875"
Front Wheel/Rear Wheel	Front: 16"   Rear: 12"
Weight	35 lbs.
Weight Limit	350 lbs
Warranty	Limited Lifetime