

## Aluminum Crutches with Accessories



- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lbs on aluminum crutch (excluding pediatric)

**PRODUCTS****Description**

10400-8	Adult, 8 pr/cs
10401-8	Youth, 8 pr/cs
10402-8	Tall Adult, 8 pr/cs
10416-1	Pediatric, 1 pr/cs

**SPECIFICATIONS****Description**

Weight Limit (10400   10401   10402)	350 lbs.
Weight Limit (10416)	175 lbs.
Warranty	Limited Lifetime