



Deluxe Bariatric Shower Chair with Cross-Frame Brace With Tool-free Removable Back

- Frame has width/depth cross braces attached with aircraft-type rivets
- Large seat with 16.5" depth
- Seat height adjusts in 1" increments
- Durable, blow-molded plastic bench
- Tool-free removable back
- Aluminum frame is lightweight, durable and corrosion proof
- Drainage holes reduce slipping
- Weight Capacity 600 lbs

PRODUCTS

Description

12029-2	Assembled, 2/cs
---------	-----------------

SPECIFICATIONS

Description

Outside Legs	21.25" (W) x 18.5" (D)
Seat Dimensions	16.5" (W) x 16.5" (D)
Seat Height	16.5" - 21.5"
Weight	9 lbs.
Carton Shipping Weight (-1 -2)	10 20 lbs.
Weight Limit	600 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12028T	Replacement Tips for 12028, 12029, 4/bx
--------	---