



Teak Adjustable-Height Shower Bench

Tool-free assembly

Item #	Product Description	UOM
RTL12351KDR	Rectangular Seat, Knocked Down	1/cs

- Attractive high-quality teak is renowned for its durability and weathering properties.
- Teak will naturally take on a handsome silver gray patina after several months of exposure to water.
- Ideal for use in shower stall.
- Ideal as a woman's grooming stool.
- Legs are lightweight, durable, and height adjustable in 1" increments.
- Easy, tool-free assembly.
- Large rubber tips reduce slipping.

specifications

Outside Legs	16.25" (W) x 16.25" (D)
Seat Dimensions	19.75" (W) x 15.25" (D)
Seat Height	13" - 18"
Weight	6 lbs.
Carton Shipping Weight	9 lbs.
Weight Limit	300 lbs.
Warranty	Limited One Year

Deluxe Handheld Shower Spray with Diverter Valve

Item #	Color	UOM
12037	White	1/bx

- On/Off switch built into handle for easy access.
- Wall holder and extra-long 80" white reinforced nylon hose provide extra convenience.
- Comes with diverter valve (Figure A) that allows either the regular shower head to operate or be used as a handheld shower spray.
- Useful when shower chair, bench or stool is being used.

specifications

Warranty	Limited Lifetime
----------	------------------