

BLUE STREAK**Blue Streak**

Single Axle (K0001)

FEATURES:

- Blue Powder coated frame, black cross brace.
- Flip back, padded desk length arm.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Solid rubber tires mounted on composite wheel provides durability and low maintenance.
- Comes with push to lock wheel locks.
- Comes with calf strap.
- Comes with swing-away footrests or elevating legrests.
- Weight Capacity: 250 lbs.
- 3 Year Limited Warranty.

1-Pack Available for All Wheelchair Models

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*

Blue Streak – SINGLE AXLE

BLS18FBD-SF	18", Blue Streak, Flip Back Desk Arm, Swing-away Footrest	1/cs
BLS18FBD-ELR	18", Blue Streak, Flip Back Desk Arm, Elevating Legrest	1/cs
BLS20FBD-SF	20", Blue Streak, Flip Back Desk Arm, Swing-away Footrest	1/cs
BLS20FBD-ELR	20", Blue Streak, Flip Back Desk Arm, Elevating Legrest	1/cs

REPLACEMENT FRONT RIGGINGS (All Silver Sport Wheelchairs Include Front Riggings.)

STDS5N22SF	Swing-Away Footrest	1pr/cs
BS-ELR	Swing-Away Elevating Legrest	1pr/cs

ARMS

STDS5N22DR	Right, Flip Back Desk Arm	1/cs
STDS5N22DL	Left, Flip Back Desk Arm	1/cs

WHEEL LOCKS (Push to Lock)

STDS5N221R	Right Wheel Lock, Push to Lock, Detachable Arm	1/ea
STDS5N221L	Left Wheel Lock, Push to Lock, Detachable Arm	1/ea

ANTI TIPPERS WITH WHEELS

STDS833		
STDS833N		

SEAT BELTS

STDS850 (Auto Clasp Type) 1/1bx		
STDS851 (Velcro® Type Closure)		1/bx

BARIATRIC SEAT BELTS

STDS855 (Auto Clasp Type)		1/bx
STDS856 (Velcro® Type Closure)		1/bx

WHEEL LOCK EXTENSIONS

STDS801		1/ea
---------	--	------

ANTI FOLD BAR

STDS806		1/cs
---------	--	------

WHEELCHAIR CANE/CRUTCH HOLDER

STDS1034		1/bx
----------	--	------

HEEL LOOPS

STDS831		1pr/bx
---------	--	--------

LIMB SUPPORT

WASR, Right		1/cs
WASL, Left		1/cs

ACCESSORIES ON PAGES 60-81

1-Pack Available for All Wheelchair Models

SPECIFICATIONS

Seat Width	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (with riggings)	Product Weight	Carton Shipping Weight	Weight Limit
18"	24"	12.5"	36"	19.5"	17.75"	18"	16"	42"	41.25 lbs.	47 lbs.	250 lbs.
20"	26"	12.5"	36"	19.5"	17.75"	20"	16"	42"	41.25 lbs.	47 lbs.	250 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

To order call toll free: 877.224.0946

WHEELCHAIRS

drive™