

## Cane Seat



FIGURE A



10365



RETAIL BOX



10365

- Provides a comfortable seat to rest on when open, and a sturdy support cane when closed (Figure A)
- Manufactured with sturdy, extruded aluminum tubing
- Tripod design with vinyl contoured tipped legs
- 9" diameter seat
- 33" tall when closed (Figure B)
- Weighs 2.5 lbs
- Weight Capacity 250 lbs

### PRODUCTS

PRODUCTS	Description
10365	Cane Seat Folding, 1/cs
RTL10365-ADJ	Adjustable Height Cane Seat, 1/cs

### SPECIFICATIONS

SPECIFICATIONS	Description
Seat To Floor Height	21"
Cane Height	34" - 38"
Seat Height	19" - 22"
Weight Limit	250 lbs.
Warranty	Limited Lifetime