

Cane/Sling Seat



FIGURE A



FIGURE C

- Provides a comfortable seat to rest on when open, and a sturdy quad cane when closed (Figure A)
- Manufactured with sturdy, extruded aluminum tubing
- Four sturdy legs with vinyl contoured tips
- Nylon sling bench, 15" x 5.75"
- Foam grip handles
- 34" tall when closed
- Seat to floor height: 20.5"
- Weighs 2.5 lbs
- 250 lb Weight Capacity

PRODUCTS

Description

RTL10360

1/cs

SPECIFICATIONS

Description

Dimensions

15" (W) x 5.75" (H)

Height Closed

34"

Seat to Floor Height

20.5"

Weight Capacity

250 lbs.

Warranty

Limited Lifetime on frame, 6 Month Warranty on Sling