



Bariatric Quad Canes

Offset Aluminum Canes with Tab-Loc Silencer

Item #	Product Description	Color	UOM
10316-4	Small Base, Vinyl Contoured Hand Grip	Chrome	4/cs
10317-4	Large Base, Vinyl Contoured Hand Grip	Chrome	4/cs

Item #	Color	UOM
10381SL-6	Silver	6/cs
10381BLK-6	Black	6/cs

- 500 lb. Weight Capacity.
- Lightweight, extra-strong base.
- Manufactured with aluminum tubing, steel base.
- Comfortable vinyl, contoured hand grips.
- Four point base, combined with offset handle provides additional stability and support.
- Aluminum handle, reinforced steel base. (Figure A)
- Handle height adjusts from 30" to 39".

- Tab-Loc silencer provides added safety to ensure a locked height adjustment. (Figure A)
- Features a padded grip. The flat side of the grip faces up to reduce stress that the user's hand may experience.
- Ergonomically designed padded handle.
- Comes standard with wrist strap.
- Manufactured with sturdy, extruded aluminum tubing.
- Handle height adjusts from 28.75" to 37.75".
- 300 lb. Weight Capacity.

specifications

Handle Height	30" - 39"
Weight Limit	500 lbs.
Warranty	Limited Lifetime

specifications

Handle Height	28.75" - 37.75"
Weight Limit	300 lbs.
Warranty	Limited Lifetime