

Item #	Seat	UOM
CPN22FBA	22" Sling Seat	1/cs



HCPCS
K0824



Cirrus Plus HD
Heavy-Duty Folding Power Wheelchair
Rear-Wheel Drive

- Foldable, lightweight frame for easy portability. (Figure A)
- Steel double cross-brace frame.
- Adjustable-tension back upholstery.
- Armrests are padded, adjustable and removable.
- Adjustable-length controller mount.
- Adjustable-tension back upholstery.
- Swing-away footrests with heel loops.
- Optional elevating leg rests (Item #CPELR) sold separately.
- Optional 2" seat depth extension (Model #STDS1818N) sold separately.
- Positioning belt included.
- Padded back upholstery with storage pouch and cushioned seat.
- Programmable controller.
- Flat-free, non-marking tires.
- Easy freewheel operation.
- Manual wheel locks.
- General foam seat cushion included.
- Anti-tip wheels.

specifications

HCPCS	K0824	Batteries (pair)	12V x 36AH (U1)
Weight Capacity	400 lbs.	Charger	4A Offboard
Top Speed	5 mph	Brakes	Electromagnetic
Maximum Range	15 miles	Caster Wheels	8" x 2" Solid
Turning Radius	35"	Drive Wheels	12.5" x 2.25" Flat-free
Climbing Angle	6°	Base Weight	90 lbs.
Ground Clearance	4"	Battery Weight (pair)	56 lbs.
Dimensions	29.5" (W) x 42.5" (L)	Seat Weight	N/A
Floor-to-Seat Height	20"	Total Weight	146 lbs.
Seat Dimensions	22" (W) x 16" (D)	Warranty on Frame	Lifetime
Controller	PG VSi 50A	Warranty on Electronic Controller/Drive Train Components	14 Months
Motors (2)	24V x 320W x 4600 rpm	Warranty on Batteries	6 Months