


Clamp-On Tub Rail



FIGURE A



RETAIL BOX



12036



12041

- Durable, white powder-coated steel construction
- 14.5" - 17" height adjustment (for RTL12036-ADJ only)
- Width adjusts from 3" - 7"
- Stainless steel locking mechanism and non-wearing rubber pads provide scratch-proof security
- Fits most tubs (not for use with fiberglass tubs)
- Tool-free installation

PRODUCTS

Description

RTL12036-ADJ	Adjustable Height, Parallel to the Tub, 1/bx
12036	Non-Adjustable Height, Parallel to the Tub, 1/bx

SPECIFICATIONS

Description

Adjustable Clamp	3" (smallest) 7" (widest)
Weight (12036 12041)	6.5 lbs.
Weight (RTL12036-ADJ)	7.0 lbs.
Height (RTL12036-ADJ)	14.5" - 17"
Height (12036 12041)	14.5"
Carton Shipping Weight (12036 12041)	8.5 lbs.
Carton Shipping Weight (RTL12036-ADJ)	9 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime