



Deluxe Folding Travel Walker with 5" Wheels and Fold-Up Legs



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



RETAIL BOX

- Easily folds tool-free to ½ the size of a traditional walker (Figure A - B) and can be stored in a "Ready Set Go" Pouch (Figure C)
- Legs are attached with strong elastic cords so legs can be folded (Figure D)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces
- Two side-storage pouches with large pockets for personal items are held in place by hook and loop fastener tabs to both side frames

PRODUCTS

PRODUCTS	Description
RTL10263KDR	Adult, Flame Red Crossbar with Red Frame, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Width	17" inside hand grip 22" inside back legs
Overall Width	26.5"
Depth (Opened at Base)	19.75"
Height	29.5 - 37"
Weight	8 lbs.
Carton Shipping Weight	10 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

PARTS

PARTS	Description
10263DRPOUCHES	Side Pouches
10263RL	Rear Leg & Glide
10263DKRWHEEL	Front Leg & Wheel
10263LEGCLIP	Leg Clips
10263DRCARRYBAG	Carry Bag