



## Deluxe, Tilt-Top Overbed Table



FIGURE A



FIGURE B

- Easily assembled
- Walnut, wood grain low-pressure laminate top
- Top can be tilted 33° in either direction (Figure A)
- Top can be raised or lowered in infinite settings between 29.5" - 46"
- Swivel casters allow for easy maneuverability
- Chrome-plated steel "H" base provides security and stability
- Tabletop locks securely when height adjustment handle is released It can be raised with the slightest upward pressure
- For proper assembly, do not remove bottom screw (Figure B)

PRODUCTS	Description
13008	1/cs

SPECIFICATIONS	Description
Height Adjustable	29.5" - 46"
Weight	21 lbs.
Carton Shipping Weight	23 lbs.
Weight Limit	50 lbs.
Base Dimensions	15.5" (D) x 26.5" (W)
Tabletop Surface Dimensions	15" (D) x 30" (W)
Warranty	Limited Lifetime

ACCESSORIES	Description
13003NC	1/ea
13003NCL	1/ea

PARTS	Description
13008H	Hardware, 1/ea
13007CL	Caster Locking, 1/ea
13007C	Caster Non-Locking, 1/ea
13008AM	Adjusting Mast, 1/ea
13008B	Plastic Bushing, 1/ea
13008K	Adjustment Knob, 1/ea
13008AB	Adjustment Bolt Bushing, 1/ea
13008EL	Elevating Lever, 1/ea
13008AAB	Adjustment Bolt, 1/ea
13008-T	Table Top, 1/ea
13008TB	Top Mounting Bracket, 1/ea
13008CP	Caster Plug (4 ea), 1/bx