

Durable 4 Wheel Rollator with 7.5" Casters

HCPCS: E0143+E0156



- Durable steel frame
- Convenient mobility with a comfortable, restful seat
- Pouch under seat for added privacy and security of personal items
- Easy-turn lever for quick height adjustment
- Large 7.5 inch caster wheels are great for indoor and outdoor use
- Deluxe hand brake can be pressed when needed or pushed down to lock rollator into place
- Removable, hinged, padded backrest can be folded up or down
- Easily folds for transport or storage

PRODUCTS	Description
10257BL-1	Rollator,Steel,4 Wheel,Blue,1/ea
10257RD-1	Rollator,Steel,4 Wheel,Red,1/ea

SPECIFICATIONS	Description
Carton Shipping Weight	22.2 lbs.
Casters	7.5"
Handle Height	31" - 35"
Inside Hand Grip Width	17"
Number of Wheels	4
Seat Depth	12"
Seat Width	12"
Seat to Floor Arm Height	21"
Unit Dimensions	23.5" (W) x 25.5" (L)
Warranty	Limited Lifetime
Weight	18.6 lbs.
Weight Capacity	300 lbs

PARTS	Description
9502F1025714	Right Handlebar, 10257, 1/ea, EA
9502F1025715	Left Handlebar, 10257, 1/ea, EA
95032F260N06	Dust Cap,260N,10257, 1/pr, PR
9502F1025712	Padded Seat W/Brackets, 10257, 1/ea, EA
9502F1025702	SeatBracketsW/Screws,HardSeat,10257,1/pr, PR
9502F1025713	SeatBrackets&Hrdwr,PaddedSeat,10257,1/pk, PKG
9502F1025711	WheelBrakeRight10257,1ea, EA
9502F1025710	WheelBrakeLeft10257,1ea, EA

PARTS

Description

9502F1025707	FoldingHingeAssembly10257,1ea, EA
9502F1025709	FrontAxleBolt10257,1ea, EA
95032F260N10	Height Adjustment Knobs,260N, 10257 1pr, PR
95032F260N09	Wheel W/Bearings, Grey, 260N & 10257,1ea, EA
95032F260N09-BLK	2F Wheel W/Bearings, 260N/10257,1ea, EA
9502F1025708	Fork10257,Blue,1ea, EA
9502F1025708R	RedFork10257,1ea, EA
9500004101	Cable,Fastner,1/ea, EA
9501S72603	Backrest Pivot Brackets for R726,1 pr, PR
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS001000	Axle Cap 8mm, 1ea, EA