

Folding Shower Chair



FIGURE A



FIGURE B



RETAIL BOX

- Product available with back (main image) and without (Figure A)
- Conveniently folds flat (Figure B)
- Blow-molded bench provides comfort and strength
- Drainage holes in seat reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Angled legs provide additional stability
- Retail packaged

PRODUCTS

Description

12486	Without Back, 1/cs
-------	--------------------

SPECIFICATIONS

Description

Outside Legs (12486)	20.5" (W) x 19.25" (D)
Outside Legs (12487)	17.75" (W) x 14" (D)
Seat Dimensions (12486)	19.75" (W) x 12" (D)
Seat Dimensions (12487)	18" (W) x 11" (D)
Seat Height (12486 12487)	17.5" - 18.5" 16.25" - 18.25"
Weight (12486 12487)	4.5 lbs. 5 lbs.
Carton Shipping Weight (12486 12487)	6.6 lbs. 7 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12486T	Replacement Tips for 12486, 12487, 4/bx
--------	---