



Shown with Half Rails.
Sold Separately

accessories

15030HBL Sold Separately



A

Full-Electric Bariatric Bed, 42"

Item #	Product Description	UOM
15300	Bariatric Bed	1/cs
15300BV-1HR	Bariatric Bed with 1 pair "T" Rails (item #15021BV)	1/cs
15300BV-PKG	Bariatric Bed Package (contains Bariatric Bed, 42" Bariatric Foam Mattress, item #15301, and 1 pair "T" Rails, item #15021BV)	1/cs

- Quiet, smooth operation.
- Heavy-duty frame ensures strength and patient safety.
- All steel construction.
- Larger sleep surface than conventional bed.
- Split-pan design with removable bed ends is easy to set up.
- No tools necessary for set up.
- Hand control (included) provides for multiple bed positioning for patients.
- Low air loss/alternating pressure mattress available.
- "T" rails available. (Sold separately) (Figure A)
- Emergency manual crank included.
- Weight Capacity 600 lbs.

specifications

Bed Height	18.25" - 26.25"
Unit Dimensions	42" (W) x 88" (L)
Sleeping Surface	80"
Weight	231 lbs.
Carton Shipping Weight	276 lbs.
Weight Limit	600 lbs.
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components

accessories

15301	80" (L) x 42" (W) x 6" (H)	1/cs
-------	----------------------------	------