

## Standard Telescoping Full-Length Side Rail, Brown-Vein Finish




FIGURE A



FIGURE B




FIGURE C

- Attractive brown-vein finish (Figure A - B)
- Constructed of 1" steel
- Telescopes with movement of bed
- Easily adjusts up or down with spring-loaded release
- Spring-loaded crossbars are easy-to-install and remove (Figure C)

### PRODUCTS

#### Description

15001ABV	Bed Rails Full, 1pr/cs
----------	------------------------

### SPECIFICATIONS

#### Description

Warranty	Limited Lifetime
----------	------------------

### ACCESSORIES

#### Description

15050BV	Universal Crossbars, 1 pr/cs
---------	------------------------------