



3" Gel/Foam Gel-U-Seat™ Cushion

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
8046-1	3" Gel/Foam Cushion (Replaces Mason # 8045-1)	16" (W) x 16" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-2	3" Gel/Foam Cushion (Replaces Mason # 8045-2)	18" (W) x 16" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-3	3" Gel/Foam Cushion (Replaces Mason # 8045-3)	20" (W) x 16" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-4	3" Gel/Foam Cushion	20" (W) x 18" (D) x 3" (H)	E2603	1/ea, 4/cs
8046-5	3" Gel/Foam Cushion	22" (W) x 18" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-6	3" Gel/Foam Cushion	22" (W) x 20" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-7	3" Gel/Foam Cushion	24" (W) x 18" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-8	3" Gel/Foam Cushion	24" (W) x 20" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-9	3" Gel/Foam Cushion	26" (W) x 18" (D) x 3" (H)	E2604	1/ea, 4/cs
8046-10	3" Gel/Foam Cushion	26" (W) x 20" (D) x 3" (H)	E2604	1/ea, 4/cs

- The Gel-U-Seat™ cushion is specifically designed to redistribute pressure to assist in the prevention of pressure ulcers while providing patient comfort.
- Cushion incorporates high-density foam encapsulating a dual-chamber gel bladder.
- Water-resistant, low-shear stretch cover with tie straps to secure to wheelchair. The vinyl bottom is waterproof.

specifications

Cover Material	Fluid Resistant Stretch Top Waterproof Vinyl Bottom
Weight Capacity	Up to 20" 275 lbs., 22" and up 500 lbs.
Warranty	1 Year Limited