



General Use Gel/Foam Gel "E" 2" Seat Cushion

Gel-U-Seat™ Lite General Use 2" Gel Cushion with Stretch Cover

Item #	Product Dimensions	Product HCPCS	UOM
14888	18" (W) x 16" (D) x 2" (H)	E2601	1/ea

- The Gel "E" is a cushion every wheelchair user can use to redistribute pressure which helps prevent, treat and manage pressure ulcers.
- High-Density foam encapsulates a viscous gel bladder that provides proper immersion and envelopment around the bony prominences.
- The zippered, removable and replaceable cover is made of a fluid-proof, urethane-coated nylon top cover and fluid-proof vinyl bottom cover for durability.

specifications

Cover Material	Fluid-Resistant Nylon Waterproof Vinyl Bottom
Weight Capacity	250 lbs.
Warranty	1 Year Limited

Item #	Product Dimensions	Product HCPCS	UOM
8040-1	16" (W) x 16" (D) x 2" (H)	E2601	1/ea, 4/cs
8040-2	18" (W) x 16" (D) x 2" (H)	E2601	1/ea, 4/cs
8040-3	20" (W) x 16" (D) x 2" (H)	E2601	1/ea, 4/cs
8040-4	20" (W) x 18" (D) x 2" (H)	E2601	1/ea, 4/cs
8040-5	18" (W) x 18" (D) x 2" (H)	E2601	1/ea, 4/cs
8040-6	22" (W) x 16" (D) x 2" (H)	E2602	1/ea, 4/cs
8040-7	22" (W) x 18" (D) x 2" (H)	E2602	1/ea, 4/cs
8040-8	24" (W) x 16" (D) x 2" (H)	E2602	1/ea, 4/cs
8040-9	24" (W) x 18" (D) x 2" (H)	E2602	1/ea, 4/cs
8040-10	20" (W) x 20" (D) x 2" (H)	E2602	1/ea, 4/cs

- The Gel-U-Seat™ Lite cushion is specifically designed to redistribute pressure to assist in the prevention of pressure ulcers while providing patient comfort.
- Cushion incorporates high-density foam encapsulating a dual-chamber gel bladder.
- Fluid-resistant, low-shear stretch cover with tie straps to secure to wheelchair. The vinyl bottom is waterproof.

specifications

Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom
Weight Capacity	Up to 20" 275 lbs., 22" and up 500 lbs
Warranty	1 Year Limited