

Lumbar Support General Use Wheelchair Back Cushion



HCPCS: E2611



INSET B



14906

- Memory-foam back cushion provides support, comfort and a tactile curve to improve posture
- Durable Masonite board prevents hammocking in sling-back wheelchair Easily removed as per patient preference
- Conforms to lumbar curve to support lower back
- Sealed liner protects foam core from moisture
- Non-slip vinyl back prevents the cushion from sliding and is waterproof
- Adjustable Straps easily attach to wheelchair
- Meets flammability standard TB117

PRODUCTS

	Description
14906	Wheelchair Back Cushion, E2611, 1/ea
14889	Wheelchair Back Cushion, E2611, 1/ea
14920	Wheelchair Back Cushion, E2611, 1/ea
8033	Wheelchair Back Cushion, E2612, 1/ea

SPECIFICATIONS

	Description
Top Fabric Material	Black PU (stretch fabric)
Bottom Fabric	PVC + Knitted Fabric
Warranty	18 Months Limited