





**Go-Lite Bariatric Steel Rollator,
Padded Seat, Loop Locks**

Item #	Color	UOM
10215BL-1	Blue 	1/cs
10215RD-1	Red 	1/cs

- 8" casters with soft-grip tires are ideal for indoor and outdoor use.
- Soft padded oversized seat.
- Comes with large 21" (W) x 6" (H) x 9" (D) basket that can be mounted under seat.
- Tool-free removable padded backrest for comfort.
- Special loop lock made of internal aluminum casting operates easily and ensures safety. (Figure A)
- Handles easily adjustable in height with self-threading knob. (Figure B)

specifications	
Handle Height	35.25" - 39.5"
Unit Dimensions	30.5" (W) x 26.75" (L)
Seat Dimensions	18" (W) x 22" (H) x 13" (D)
Weight	26 lbs.
Carton Shipping Weight	30.5 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime