



“T” Style Half Rails
with Brown-Vein Finish

Half-Length Siderail
No-Gap Style with Brown-Vein Finish

Item #	Color	UOM
15021BV	Brown-Vein	1 pr/cs

Item #	Color	UOM
15201BV	Brown-Vein	1pr/cs

- Ideal protection without confining the patient.
- Provides maximum protection when bed is elevated.
- Rail easily attaches to bed. (Figure A - B)
- For use with “Pan-style” beds and Bariatric Beds.
- Attractive brown-vein finish that is easy to maintain and less institutional looking. (Figure C)
- Constructed of 1" steel.
- For use with most manufacturers' beds.

- Unique design prevents individuals from being lodged between bars.
- Ideal for maximum protection when bed is elevated.
- Fits most hospital beds with steel frame. (Figure A)
- Attractive brown-vein finish is easy to maintain and less institutional looking. (Figure B)
- Constructed of 1" steel.
- Not for use with “Pan-Style” bed surface.
- Easy to install.

specifications

Warranty Limited Lifetime

specifications

Warranty Limited Lifetime

accessories

14198	Clamp for 15201BV Rails (come standard)	1 pr/cs
14195	Clamp for 15202 Rails (come standard)	1 pr/cs