

Hinged Toilet Seat Riser

HCPCS: E0244



INSET A

INSET B

- Raises the height of the toilet 3.5" for individuals that require assistance sitting on and standing from the toilet
- Hinge design allows riser to be moved out of the way when not in use (Figure A)
- Heavy-duty molded plastic construction provides strength and durability
- Lightweight, portable, and easily attaches to toilet bowl
- Easy-to-install on most standard toilets
- Two models to accommodate either standard or elongated (Figure B) toilet bowls

PRODUCTS

Description

RTL12607	Standard Seat, 1/ea
RTL12608	Elongated Seat, 1/ea

SPECIFICATIONS

Description

Seat Dimensions (12607)	13.5" (W) x 17.5" (D) x 3.5" (H)
Seat Dimensions (12608)	14" (W) x 19.25" (D) x 3.5" (H)
Weight (12607 12608)	5.07 lbs. 5.62 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime