

## Lightweight Bariatric Homecare Bed

HCPCS: 15300LW: E0301 | 15302LW and 15303LW: E0302



- Lightweight frame design offers providers easier and more manageable bed deliveries and assemblies
- Allows for additional space for patient comfort and recovery as well as caregiver efficiency
- Product Safety & Performance Standards: Meets IEC 60601-2-52; 60601-1; 60601-1-2; 60601-1-11
- Safe working load capacity of 600 lbs
- Full electric bariatric homecare bed with dual motor assembly for head and frame adjustability
- Ultra lightweight frame is easily set up by 1 individual; no links to connect and no tools required
- Square tube steel frame slat deck provides durable strength, no-sag mattress support, and is easy to clean
- Scratch resistant wood grain panel bed ends include hi lo motors for frame height adjustability
- UL approved motor is self contained for reduced weight and noise and can be easily installed and removed with patient in bed
- Reduce service calls with large glow in the dark hand pendant that is easy-to-use for adjustments; LED indicators available for power and lockout included
- Battery Back up accessory available (#15302LWB) and sold separately
- Half Rails accessory available (15021BV-LW) and sold separately
- Rail accessories meet FDA and ASTM entrapment guidelines and standards
- 
- 

### PRODUCTS

PRODUCTS	Description
15300LW	Lightweight Full Elec 42" Bariatric Bed
15302LW	Lightweight Full Elec 48" Bariatric Bed
15303LW	Lightweight Full Elec 54" Bariatric Bed

### SPECIFICATIONS

SPECIFICATIONS	Description
Battery	Backup available but not included
Bed Height	15300LW: 34" - 42" with casters   15302LW and 15303LW: 35" - 44" with casters
Carton Shipping Weight	15300LW: 255 lbs (2 boxes)   15302LW: 294 lbs (3 boxes)   15303LW: 312 lbs (3 boxes)
Casters	5" casters (2 locking and 2 unlocking)
Color	Brown
Ground Clearance	15300LW: 18 - 26"   15302LW and 15303LW: 16 - 24"
Motor/Drive System	15300LW: Dual Full Electric w/Hi Lo   15302LW and 15303LW: Linear Full Electric w/Hi Lo
Sleep Surface	15300LW: 42"W x 80" L   15302LW: 48"W x 80" L   15303LW: 54"W x 80" L
Warranty	5 years on Motor and Hand Control   Limited Lifetime on Welds and Frame   1 year on all other parts and components
Weight Capacity	15300LW: Safe working load 600lbs (to include occupant and any accessories)   15302LW: Safe working load 750lbs (to include occupant and any accessories)   15303LW: Safe working load 1000lbs (to include occupant and any accessories)

**SPECIFICATIONS****Description**

Material Composition

Bed Frame: Steel | Bed Ends: Steel base &amp; Wood grain panels | Motors: electrical components

**ACCESSORIES****Description**

15302LWB

LW Bari Bed Accessory battery pk,48/54, BX

15021BV-LW

Half Side Rail,42in/48/54 LW Baribed,1pr, BX