



Travelite Transport Chair **Super Light, Folding Transport Chair with Carry Bag and Flip-Back Arms**

Item #	Seat Width	UOM
TC005GY	18"	1/cs

- The Travelite Chair in a Bag is a lightweight yet sturdy aluminum transport chair in gunmetal gray.
- It is easily transportable and comes with a compact carry bag with wheels making it ideal for traveling. (Figure A - B)
- Folds for storage and transportation. (Figure C)
- It is ideal for use indoors and outdoors.
- 18" square nylon seat comes with lift-up straps.
- Comes with wheel locks. 8" solid tires.
- Comes with padded, flip-up desk arms and height adjustable footrests.

specifications	
Overall Dimensions	22.5" (W) (open) 11.5" (W) (closed) x 37.5" (H)
Seat-To-Floor Height	19"
Back Height	17.5"
Seat Dimensions	18" (W) x 18" (D)
Overall Length (with riggings)	33"
Weight	20.5 lbs.
Carton Shipping Weight	23.2 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime

Item #	Color	Seat Width	UOM
SL18	Silver Frame and Blue Upholstery	18"	1/cs

- Padded flip-back armrests. Fold-down footrests.
- Weighs only 19 lbs. Easy to push or transport.
- Aluminum frame with lightweight nylon upholstery.
- Composite wheels are lightweight and maintenance free.
- Folds like a conventional folding chair. (Figure A)
- Comes with carry bag (item #835) with shoulder strap. (Figure B)
- Standard rear wheel locks and seat belt.

specifications	
Overall Dimensions	23" (W) (open) x 37" (H)
Seat-To-Floor Height	20"
Back Height	17"
Seat Dimensions	18" (W) x 17" (D)
Overall Length	30"
Weight	19 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime

accessories		
835N	Replacement Carry Bag	1/ea