

Nitro Duet Rollator and Transport Chair

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F



FIGURE G



FIGURE H

- Combines the features of a rollator and transport chair into one unit. The user can ambulate independently or be pushed safely by a caregiver
- The tool-free two-position contoured backrest is reversible depending on whether you are using the Nitro Duet as a rollator or transport chair
- Back support height easily adjusts with tool-free buttons
- Caster fork design and large 10" front casters enhances turning radius and rolling comfort
- Cross-brace design allows for side-to-side folding and added stability
- Brake cables routed inside the frame add safety and convenience

PRODUCTS

Description

| | |
|------------|------------------------|
| RTL10266DT | Nitro Duet, Red, 1 Ea. |
|------------|------------------------|

SPECIFICATIONS

Description

| | |
|------------------------|------------------|
| Front Wheels | 10" |
| Handle Height | 33" 37" |
| Inside Hand Grip Width | 18" |
| Number of Wheels | 4 |
| Rear Wheel | 8" |
| Seat Depth | 14" |
| Seat To Floor Height | 20.5" |
| Seat Width | 18" |
| Warranty | Lifetime Limited |
| Weight Capacity | 300 lbs. |