



Nitro Elite CF, Carbon Fiber Rollator

Item #	Product Description	Color	UOM
RTL10266CF	Nitro Elite CF	Black Carbon Fiber	1/ea

- Luxurious super lightweight carbon fiber frame. (Figure A)
- Weighs only 12 lbs.* for powerful performance with swift and effortless mobility.
- Opulent and removable zippered storage bag. Unique bag attachment keeps bag securely in place when rollator is open or folded. (Figure B)
- Brake cable inside frame for added safety. (Figure C)
- Handle height easily adjusts with unique push button.
- Caster fork design enhances turning radius.
- Large 10" front casters allow optimal steering and rolling comfort. (Figure D)
- Easily fold with one hand to ultra-compact size for storage. Handle is located on seat. (Figure E)
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable.
- EZ clip cane holder with HurryCane™ compatible platform.

specifications

Handle Height	33.25" - 38.2"
Unit Dimensions	23" (W) x 27.75" (L)
Seat Dimensions	17.75" (W) x 20.75" (H) x 8.5" (D)
Weight	12 lbs.*
Carton Shipping Weight	19 lbs.
Weight Limit	300 lbs.
Warranty Limited Lifetime	Frame
5 Year Limited Warranty	Brake and Handle Housing

*excluding bag

accessory

10266-CH	Nitro Cup Holder	1/ea
----------	------------------	------