

## Nitro 3-Wheel Rollator

**PreserveTech™**  
Stays Cleaner • Lasts Longer



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F



FIGURE G




FIGURE H

The Nitro 3-Wheel Rollator is the newest addition to our premium line of Nitro rollators. Featuring a hidden inner brake cable, the functional and modern Nitro 3 is designed with safety in mind while indoors or outdoors. The Nitro 3 easily navigates narrow spaces and folds flat while in the standing position so it can roll, park and store anywhere. Also included is a large, removable pouch to keep personal items close. The Nitro 3 comes pre-assembled, so you can unbox, unfold and get your customers going!

- Stylish and modern design features lightweight aluminum frame with built-in brake cable for extra safety and protection
- Ergonomic hand grips provide comfort and control, push down on handbrakes to engage parking brake feature
- Quickly and easily adjust the height of the handles with the push-button height adjustment
- Large, removable storage pouch is perfect for your personal belongings, or transporting items to and from
- Easy side to side fold allows the rollator to stay in the standing position, making it easy to roll, park and store anywhere
- Dual front fork design enhances turning radius while providing added safety
- Large 9" front wheel allows for optimal steering and conveniently moves over most kinds of terrain, whether indoors or out
- Out of box assembled
- Drive PreserveTech products are manufactured with an antimicrobial\* agent that helps to minimize the growth of bacteria, mold and mildew on treated surfaces. Treated surfaces on this product include the handgrips

PRODUCTS	Description
RTL10266TWHL	Nitro 3-Wheel Rollator, Red, 1/EA
SPECIFICATIONS	Description
Dimensions	22" W x 29" D
Front Wheels	9"
Handle Height Range	32.5" - 37.5"
Inside Hand Grip Width	14"
Number of Wheels	3
Product Weight	12.7
Rear Basket Size	12.5" W x 13.5" H x 11.5" D
Rear Wheels	8"
Warranty	Limited Lifetime
Weight Capacity	300 lbs
Width inside back legs	21.5"