

Nitro HD Rollator



INSET A




INSET B



INSET C



INSET D



INSET E



INSET F

- Attractive, Euro-style design
- Brake cable inside frame for added safety (A)
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort (B)
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius (C)
- Large 10" front casters allow optimal steering and rolling comfort (C)
- Easily folds with one hand to ultra-compact size for storage (D)
- Under seat locking mechanism holds the Nitro closed when frame is folded (E)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (F)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

PRODUCTS

PRODUCTS	Description
10266HD-BK	Black, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Handle Height	36" - 41"
Unit Dimensions	28" (L) x 28.25" (W)
Seat Dimensions	23" (W) x 10" (D)
Seat Height	23.5"
Weight	20.5 lbs.
Weight Limit	450 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake Handle & Housing)	5 Years

ACCESSORIES

ACCESSORIES	Description
892	Cane Holder, 1/ea
10266-CH	Nitro Cup Holder, 1/ea