



Offset Handle Aluminum Canes Foam Rubber Grip



FIGURE A



FIGURE B




FIGURE C




FIGURE D




RTL10303CF



RTL10303PF



RTL10303SB



- Ergonomically designed handle with soft Foam Grip provides comfort and security
- Comes standard with foam rubber grip and wrist strap (Figure A)
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring Locking ring prevents rattling (Figure B)
- Available in Silver (main image), Black (Figure C) or Bronze (Figure D)
- Handle height adjusts from 30" to 39"
- Weight Capacity 300 lbs

PRODUCTS	Description
10303-6	Silver, 6/cs
10303-1	Silver, 1/ea
RTL10306	Black, 1/ea
RTL10307	Bronze, 1/ea
10308C-6	Silver, Black & Bronze, 6/cs
RTL10303	Silver, 1/ea
RTL10303CF	Camouflage, 1/ea
RTL10303PF	Pink Floral, 1/ea
RTL10303SB	Stars, 1/ea
RTL10303SW	Swirl, 1/ea
SPECIFICATIONS	Description
Handle Height	30" - 39"
Weight Limit	300 lbs.
Warranty	Limited Lifetime