

SLINGS FOR FLOOR LIFTS



Divided Leg Sling




Toileting Sling




Item #	Sling Points	Chain/Strap	Description	Design	Material	Size	Product Dimensions	Wt. Cap.	Cradle Points
13262P	4	Not Required	Divided Leg Sling	Solid	Nylon	Petite	55" x 34"	450 lbs.	4 or 6
13262S	4	Not Required	Divided Leg Sling	Solid	Nylon	Small	60" x 38"	450 lbs.	4 or 6
13262M	4	Not Required	Divided Leg Sling	Solid	Nylon	Medium	62" x 39"	450 lbs.	4 or 6
13262L	4	Not Required	Divided Leg Sling	Solid	Nylon	Large	63" x 44"	450 lbs.	4 or 6
13226T	4	Not Required	Toileting Sling	Solid	Nylon	Large	33" x 35"	450 lbs.	4 or 6

SLINGS FOR STAND ASSIST LIFTS



Stand Up Sling

Transfer Sling



Item #	Sling Points	Chain/Strap	Description	Design	Material	Size	Product Dimensions	Wt. Cap.	Cradle Points
13263A	2	Not Required	Stand-Up	Solid	Nylon	Large	44" x 16"	440 lbs.	2
13263B	2	Not Required	Stand-Up	Solid	Nylon	X-Large	55" x 16"	440 lbs.	2
13264A	2	Not Required	Transfer	Solid	Nylon	Large	37" x 34"	350 lbs.	2
13264B*	2	Not Required	Transfer	Solid	Nylon	X-Large	47" x 44"	440 lbs.	2



Sit to Stand Sling



Sit to Stand Strap

Item #	Sling Points	Chain/Strap	Description	Design	Material	Size	Product Dimensions	Wt. Cap.	Cradle Points
13263D	2	Not Required	Sit to Stand Padded Sling	Solid	Nylon	Small	26" – 42"	400 lbs.	2
13263E	2	Not Required	Sit to Stand Padded Sling	Solid	Nylon	Medium	30" – 48"	400 lbs.	2
13263F	2	Not Required	Sit to Stand Padded Sling	Solid	Nylon	Large	36" – 60"	400 lbs.	2
13263G	2	Not Required	Sit to Stand Padded Sling	Solid	Nylon	X-Large	50" – 79"	600 lbs.	2
13263H	2	Not Required	Sit to Stand Buttock Strap	Solid	Nylon	-	71"	400 lbs.	2

Drive slings can be used with any competitive floor-style lift.

Sling Chains (13019), Bariatric Chains (13019B) or Sling Straps (13232) Sold Separately