




Item #	Seat	UOM
PHOENIXHD3	17.5" Folding Seat	1/cs



HCPCS
K0801



A



B

Phoenix HD 3 Heavy-Duty Travel Scooter

- Quick and easy disassembly. (Figure A)
- Includes interchangeable color panels in Red and Blue. (Figure B)
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with two-tone upholstery.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Swivel seat with fold-down backrest.
- Easy freewheel operation.

specifications

HCPCS	K0801	Batteries (pair)	12V x 20AH
Weight Capacity	350 lbs.	Charger	1.5A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Maximum Range	15 miles	Front Wheel	9"x 3" Flat-free
Turning Radius	43"	Rear Wheels	9"x 3" Flat-free
Climbing Angle	6°	Base Weight	63 lbs.
Ground Clearance	1.75"	Battery Weight (pair)	31 lbs.
Dimensions	39" (L) x 22.5" (W)	Seat Weight	21 lbs.
Floor-to-Seat Height	20" - 22"	Total Weight	115 lbs.
Drive System	24V x 350W x 3600rpm	Heaviest Piece Weight	32 lbs. (Front Section)
Seat Type	Folding	Warranty on Frame	Lifetime
Seat Dimensions	17.5" (W) x 17" (D)	Warranty on Electronic Controller/Drive Train Components	14 Months
Motor	24V x 350W x 3600 rpm	Warranty on Batteries	6 months