



## Premier One Foam Cushion

| Item #   | Product Dimensions                                 | UOM        |
|----------|--|------------|
| M8061-PB | 14" (W) x 14" (D) x 3" (H)                         | 1/ea, 4/cs |
| M8062-PB | 14" (W) x 16" (D) x 3" (H)                         | 1/ea, 4/cs |
| M8063    | 16" (W) x 16" (D) x 3" (H) (Replaces Mason # 8063) | 1/ea, 4/cs |
| M8064    | 18" (W) x 16" (D) x 3" (H) (Replaces Mason # 8064) | 1/ea, 4/cs |
| M8065    | 20" (W) x 16" (D) x 3" (H) (Replaces Mason # 8065) | 1/ea, 4/cs |
| M8070    | 20" (W) x 18" (D) x 3" (H) (Replaces Mason # 8070) | 1/ea, 4/cs |
| M8076    | 22" (W) x 16" (D) x 3" (H) (Replaces Mason # 8076) | 1/ea, 4/cs |
| M8077    | 22" (W) x 18" (D) x 3" (H) (Replaces Mason # 8077) | 1/ea, 4/cs |
| M8078    | 22" (W) x 20" (D) x 3" (H) (Replaces Mason # 8078) | 1/ea, 4/cs |
| M8079    | 24" (W) x 18" (D) x 3" (H) (Replaces Mason # 8079) | 1/ea, 4/cs |
| M8080    | 24" (W) x 20" (D) x 3" (H) (Replaces Mason # 8080) | 1/ea, 4/cs |
| M8081    | 26" (W) x 18" (D) x 3" (H) (Replaces Mason # 8081) | 1/ea, 4/cs |
| M8082    | 26" (W) x 20" (D) x 3" (H) (Replaces Mason # 8082) | 1/ea, 4/cs |

- The Premier One Foam cushion is designed to provide comfort and stability with a simple, economical design.
- Aerated foam keeps the patient cool and dry.
- The low shear water-resistant stretch cover helps maintain sensitive skin tissue. The vinyl bottom is waterproof.

| specifications  |   |
|-----------------|---|
| Cover Material  | Fluid-Resistant Stretch Top   Waterproof Vinyl Bottom |
| Weight Capacity | 250 lbs.   400 lbs.                                   |
| Warranty        | 1 Year Limited  |