

PreserveTech Raised Toilet Seat with Bidet

HCPCS:



Both Models

RTL12C005-WH

RTL12C006-WH

- Control panel activates and adjusts water pressure and cleansing modes
- Bidet provides front or female and back cleansing depending on needs
- Retractable bidet nozzle and guard for a more sanitary experience
- Nozzle cleaning mode on the control panel
- Connects to the water supply no electrical connection or batteries needed
- Adds 4" of height to toilet for safer sitting and standing Includes lid for better hygiene and aesthetic
- Contoured for safe and comfortable elevation
- Made from lightweight and durable HDPE to support up to 400 lb.
- Fits most standard and elongated toilets
- Antimicrobial* treated raised toilet seat helps minimize the growth of bacterial odors, mold, and fungi on treated surfaces.
- *Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs, or other disease organisms.
- Ambient water temperature only
- Easily connects to the water source behind the toilet
- Choose between warm and ambient water temperatures
- Easily connects to the water source behind the toilet and under the sink
- Comes with an inlet cap to plug warm water connection for ambient water use only to accommodate bathroom layouts

PRODUCTS

Description

RTL12C005-WH	RTS w. Bidet, Ambient Water, 1/CV
RTL12C006-WH	RTS w. Bidet, Warm Water, 1/CV

SPECIFICATIONS

Description

Weight Capacity	400 Lbs
Color	White
HCPCS	E0244
Material	HDPE
Warranty	One Year Limited