
## pressure prevention







## Med-Aire Alternating Pressure Pump and Pad System

			Product	
ltem #	Product Description	<b>Product Dimensions</b>	HCPCS	UOM
14001E	Variable Pressure Pump and Pad	78" (L) × 36" (W) × 2.5" (H)	E0181	1/ea
14001EFD	Variable Pressure Pump and Pad with End Flaps (Replaces Mason # 2000PP-ADJ)	123" (L) × 36" (W) × 2.5" (H)	E0181	1/ea
14001E-220	Variable Pressure Pump (220V) and Pad	78" (L) x 36" (W) x 2.5" (H)	E0181	1/ea
14002E	Fixed Pressure Pump and Pad (Replaces Mason # 2000PP-II)	78" (L) × 36" (W) × 2.5" (H)	E0181	1/ea
14002E-220	Fixed Pressure Pump (220V) and Pad	78" (L) x 36" (W) x 2.5" (H)	E0181	1/ea

- 130 individual bubble cells for maximum comfort and effectiveness.
- Quiet pump technology alternately inflates and deflates the air cells.
- 4 LPM (Liters per Minute) pump produces consistent air flow and pressure over a 5 minute cycle time.
- Variable pressure setting on deluxe pump allows comfort setting for maximum comfort and compliance.

specifications				
Pump Power	110 VAC, 60 Hz			
Pump Airflow	4 LPM			
Pump Cycle Time	5 Minutes			
Pump Type	Compressor			
Weight Capacity	300 lbs.			
Warranty	1 Year Pump, 90 Day Pad			

- Durable, heavy gauge, premium vinyl bubble pad
  provides superior therapeutic treatment.
- Built-in bracket for mounting easily to bed.
- Available with variable pressure pump (Figure A) or Fixed pressure pump. (Figure B)

accessories				
14003	78" (L) x 36" (W) x 2.5" (H), Standard Pad Only	1/ea		
	(Replaces Mason # 2001)			
14003-EF	123" (L) $\times$ 36" (W) $\times$ 2.5" (H), Pad with End Flaps	1/ea		
	(Replaces Mason # 2001-FL)			
14001T	Beige Tubing for 14001 and 14001E	1/ea		
14005E	Variable Pressure Pump Only for 14001E	1/ea		
	(Replaces Mason # 2002-ADJ)			
14006E	Fixed Pressure Pump Only for 14002E	1/ea		
	(Replaces Mason # 2002-II)			