



Quad Canes, Large and Small Base, with Silver Vein Finish




FIGURE A




FIGURE B

- Tab-Loc silencer provides added safety to ensure a locked height adjustment (Figure A)
- Features a padded grip
- The flat side of the grip faces up to reduce stress that the user's hand may experience
- Four point base, combined with offset handle provides additional stability and support
- Attractive silver vein finish on base (Figure B)
- Ergonomically designed padded handle
- Manufactured with an aluminum shaft and a steel base to provide the perfect balance of safety and durability
- Handle height adjusts from 29" to 38"
- Weight Capacity 300 lbs

| PRODUCTS | Description |
|----------|------------------|
| 10385-4 | Large Base, 4/cs |
| 10384-4 | Small Base, 4/cs |

| SPECIFICATIONS | Description |
|----------------|------------------|
| Handle Height | 29" - 38" |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |