



## 2-in-1 Locking Raised Toilet Seat with Tool-free Removable Arms



FIGURE A



RETAIL BOX

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms can be removed or added (Figure A)
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit
- Lightweight and portable
- Easily attaches to toilet bowl, tool-free
- Easy-to-clean
- Fits most toilets
- Model RTL12027RA is retail packaged

### PRODUCTS

PRODUCTS	Description
RTL12027RA	Elev.Toilet Seat w/ RemArms RTL, 1/cs
12027RA-4BULK	Raised Toilet Seat w/RemArms, 4/cs

### SPECIFICATIONS

SPECIFICATIONS	Description
Seat Dimensions	17" (W) x 16.5" (D) x 5" (H)
Width Between Arms	18"
Weight	5 lbs.
Carton Shipping Weight (RTL12027RA   12027RA-4BULK)	8.5   23 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime