





## Aluminum Rollator, 7.5" Casters

Fold-up and Removable Back Support, Padded Seat, Loop Locks

Item #	Color	UOM
R728BL	Blue 	1/cs
R728RD	Red 	1/cs

- Comes with new seamless padded seat (Figure A). Zippered pouch under seat for added privacy and security of personal items. (Figure B)
- 7.5" Black non-marring tires are ideal for indoor and outdoor use.
- Removable, hinged, padded backrest can be folded up or down.
- Easy-to-use deluxe loop locks. (Figure C)
- Brakes with serrated edges provide firm hold. (Figure D)
- Ergonomic handles are easy to grip and are adjustable.

### specifications

Handle Height	33" - 38"
Unit Dimensions	24" (W) x 28" (L)
Seat Dimensions	14" (W) x 23" (H) x 12" (D)
Weight	14 lbs.
Carton Shipping Weight	17 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

### accessories

9501S72606	Basket	1/ea
------------	--------	------