

Item #	Seat	Battery	Color	UOM
SFSCOUT4	16.5" Folding Seat	12AH	Red/Blue	1/ea
SFSCOUT4-EXT	16.5" Folding Seat	20AH	Red/Blue	1/ea



HCPCS
K0800



A

Scout 4 Compact Travel Scooter

- Quick and easy disassembly (Figure A)
- Color-through color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy free wheel operation
- Available with 12AH or 20AH batteries
- 20AH batteries increase scooter's maximum range

specifications			
HCPCS	K0800	Charger	2A Offboard
Weight Capacity	300 lbs.	Brakes	Electromagnetic
Top Speed	4.25 mph	Front Wheels	8" x 2" Flat-free
Maximum Range	9 miles (12AH) 15 miles (20AH)	Rear Wheels	8" x 2" Flat-free
Turning Radius	53.75"	Base Weight	58 lbs.
Climbing Angle	6°	Battery Weight (pair)	20 lbs.
Ground Clearance	2.5"	Seat Weight	16 lbs.
Dimensions	42.5" (L) x 21.5" (W)	Total Weight	94 lbs.
Floor-to-Seat Height	17" - 19"	Heaviest Piece Weight (Front Section)	39 lbs.
Seat Dimensions	16.5" (W) x 16" (D)	Warranty on Frame	Lifetime
Controller	CTE 45A	Warranty on Electronic Controller/Drive Train Components	24 Months
Motor	24V x 270W x 4700 rpm	Warranty on Batteries	12 Months
Batteries	12(V) x 12(AH)		