

**Bariatric Sentra EC Heavy-Duty, Extra-Extra-Wide Wheelchair**  
Dual Cross Brace



HCPCS  
K0007



Shown with  
optional PH-SF  
swing-away footrests

Item #	Seat Width	Armrest	Product Front Riggings	UOM
STD26ECDDA	26"	Detachable Desk Arm	Sold Separately	1/cs
STD26ECDF A	26"	Detachable Full Arm	Sold Separately	1/cs
STD28ECDDA	28"	Detachable Desk Arm	Sold Separately	1/cs
STD28ECDF A	28"	Detachable Full Arm	Sold Separately	1/cs
STD30ECDDA	30"	Detachable Desk Arm	Sold Separately	1/cs
STD30ECDF A	30"	Detachable Full Arm	Sold Separately	1/cs

- 700 lb. Weight Capacity.
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable reinforced nylon upholstery.
- Steel, spoke wheels are strong and maintenance free.
- Heavy-duty, flat-free 8" front casters are adjustable in two positions.
- Extra-heavy-duty front forks.
- Reinforced side frame, caster journals and dual cross bars.
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Extra-large, padded armrests provide comfort and support.
- Riggings have extra-large, heavy-duty, aluminum foot plates that provide maximum support, are attractive and crack proof. (Must be ordered separately)
- Comes with push-to-lock wheel locks.

### specifications

Width Open   Closed with 26" Seat	34.5"   14"	Product Weight with 20" Seat	62 lbs.
Width Open   Closed with 28" Seat	37"   14.5"	Product Weight with 22" Seat	64 lbs.
Width Open   Closed with 30" Seat	39"   15"	Product Weight with 24" Seat	66 lbs.
Height	35"	Carton Shipping Weight with 20" Seat	68 lbs.
Seat-To-Floor Height	19.5"	Carton Shipping Weight with 22" Seat	70 lbs.
Back Height	18"	Carton Shipping Weight with 24" Seat	72 lbs.
Seat Depth	20"	Weight Limit	700 lbs.
Length with Riggings	43"	Warranty	Limited Lifetime

### accessories

<b>Front Riggings</b>			
PH-SF	Aluminum Foot plate	1/pr	
PHELR	Aluminum Foot plate	1/pr	
<b>Anti-Tippers with Wheels</b>			
STDS829	Anti-Tippers With Wheels	1/pr	
<b>Wheel Lock Extension</b>			
STDS801	Wheel Lock Extension	1/ea	
<b>Bariatric Seat Belts</b>			
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea	
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea	
<b>Telescoping IV Pole Attachment</b>			
STDS820	Telescoping IV Pole Attachment	1/ea	
<b>Universal Oxygen "E" Cylinder Holder</b>			
STDS803	Universal Oxygen "E" Cylinder Holder	1/ea	
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea	
<b>Wheelchair Cane/Crutch Holder</b>			
STDS1034	Wheelchair Cane/Crutch Holder	1/ea	
<b>Heel Loops</b>			
STDS831	Heel Loops	1/pr	