

## Shower Bench with Removable Padded Arms



- Removable, padded armrests provide comfort and added security
- Blow-molded bench provides comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments

PRODUCTS	Description
12440-2	2/cs
12440-1	1/cs
12440KD-1	1/cs

SPECIFICATIONS	Description
Outside Legs	18.5" (W) x 18.25" (D)
Seat Dimensions	15.75" (W) x 13" (D)
Seat Height	16" - 21"
Weight (-1   -2)	5.5 lbs.   10.5 lbs.
Carton Shipping Weight (-1   -2)	8 lbs.   14 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
12002T	Replacement Tip 4/bx