

Shower Stool



RETAIL BOX

- Aluminum frame is lightweight, durable and corrosion proof
- Impact-resistant, composite seat
- Legs are height adjustable in 1" increments, and are crack proof and tarnish resistant
- Easy, tool-free assembly

PRODUCTS

Description

RTL12004KD	Knocked Down, White, 1/cs
------------	---------------------------

SPECIFICATIONS

Description

Outside Legs	12.5" (W) x 12.5" (D)
Seat Dimensions	13" (W) x 13" (D)
Seat Height	13.5" - 21"
Weight	3 lbs.
Carton Shipping Weight	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12002T	Replacement Tip 4/bx
--------	----------------------