



Side Walker



FIGURE A

- Designed for individuals with the use of only one hand or arm
- Lighter than a walker and more stable than a cane
- Folds easily with one hand (Figure A)
- Ideal for users who do not require a walker but need a wider base for support

PRODUCTS

Description

10240-1	Adult (Hemi-Walker), 1/ea
---------	---------------------------

SPECIFICATIONS

Description

Overall Dimensions	18.5" (W) x 16" (D)
Height	28.5" - 35.5"
Weight	4 lbs.
Carton Shipping Weight	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime