

Steel Forearm Crutches



- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability
- Ortho K Handle is hard plastic

PRODUCTS**Description**

10403	Forearm Crutch Adult, 1pr/bx
10405	Forearm Crutch, 1pr/bx

SPECIFICATIONS**Description**

Weight Limit	300 lbs.
Weight	2.5 - 2.6 lbs.
Warranty	Limited Lifetime